



Test kit, RapidSilver[™] visual

Kat. číslo produktu: 2674500

ZASTARANÁ POLOŽKA
Táto položka už nie je dostupná.

Test Kit Silver, Model RapidSilver[™]. Range 0-50 ppb. Silver As Ag Visual Comparison.

Provides fast, accurate information on silver discharges

Used by surface finishers, photographic film manufacturers and processors, and others to prevent pollution and economic loss

Kit directions also include digestion and dilution procedures needed for testing silver strike bath samples

Technické údaje

Model:	NA
Model testovacej súpravy:	RapidSilver [™]
Názov metódy:	Visual Comparison
Parameter:	Silver
Počet testov:	100
Rozsah merania:	5 - 50 ppb Ag
Smallest increment:	5 µg/L
Štandardný rozsah merania:	5 - 50 ppb Ag
Štýl kufríka:	F