



## Vacuum filter microfunnel, pk/50.

**Kat. číslo produktu:** 2831500  
**EUR Cena bez DPH:** Kontaktujte nás  
**Dátum odoslania neuvedené**

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metrical with 0.45 µm pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

---

### Technické údaje

|                  |                       |
|------------------|-----------------------|
| Materiál:        | Mixed-Cellulose Ester |
| Priemer:         | 47 mm                 |
| Veľkosť balenia: | 50                    |
| Veľkosť póru:    | 0.45 µm               |