



Filter membrane, EZ-Pak, 0.45 micron, 47 mm, white, gridded, 600/pk (millipore)

Kat. číslo produktu: 2916800
EUR Cena bez DPH: Kontaktujte nás
Dátum odoslania neuvedené

Millipore EZ Pack Membranes, 0.45 µm pore size, 47mm

Dispenser provides superior convenience and assures sterility. Unique one press on the lever and a sterile membrane is presented (horizontally and inverted) for efficient and convenient assembly. The EZ Pak membranes come as a band of membranes that easily loads into the dispenser. Packaging materials are wound into the dispenser ensuring a clear operating area and quick change of bands. The membranes are Type HA, which signifies compliance with Standard Methods criteria.

Technické údaje

Materiál:	Mixed-Cellulose Ester
Priemer:	47 mm mm
Veľkosť balenia:	600
Veľkosť póru:	0.45 µm