



Be Right™



Bottle, Storage, Wide Mouth, Write-On, 1000 mL, 6/pk

Kat. číslo produktu:	2930004
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Large dedicated labeling areas help to standardize the way bottle contents are identified. A separate area is provided for name and a larger area allows for detail of contents or safety information and includes a defined space for date information. Bottles are molded of rugged HDPE in 500 and 1000 mL sizes and are available as narrow or wide mouth with leak proof polypropylene cap.

Technické údaje

Body Material:	High Density Plastic
Kapacita:	1000 mL
Materiál:	HDPE
Package size UOM:	6pk
Veľkosť balenia:	5 /bag