



Membrane filters millipore S-pak type HA, 47 mm, 0.45 µm pore size, pk/600

Kat. číslo produktu: 2936500
EUR Cena bez DPH: Kontaktujte nás
Dátum odoslania neuvedené

For membrane filtration in microbiological analysis of aqueous samples. "HA" signifies certified compliance with Standard Methods criteria with pore size of 0.45 µm. White, cellulose and mixed-ester composition. S-Pak (filter only) is sterile. Diameter: 4

Technické údaje

Materiál:	Mixed-Cellulose Ester
Priemer:	47 mm mm
Veľkosť balenia:	600
Veľkosť póru:	0.45 µm