



Bottle, storage, autoclavable, wide mouth, 125 mL pk/12

Kat. číslo produktu:	2939103
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Polyethylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121°C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Technické údaje

Body Material:	Plastic
Kapacita:	125 mL
Materiál:	PP
Package size UOM:	12pk
Veľkosť balenia:	12 /pk