



Bottle BOD, 300 mL, disposable

Kat. číslo produktu:	2943101
EUR Cena bez DPH:	Kontaktujte nás
Dodávka do 2 týždňov	

Single 300 mL disposable BOD bottle without stopper. Plastic bottle is EPA approved for BOD measurement.

Box of 100 bottles available under article number 2943100. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Technické údaje

Body Material:	Plastic
Kapacita:	300 mL
Materiál:	PET
Veľkosť balenia:	1 each