



Large custom kit case assembly (574 x 325 x 295 mm)

Kat. číslo produktu:	4766600
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Designed to withstand the harshest conditions, this rugged, blue polypropylene case assembly includes a foam insert to prevent kit components from shifting, breaking or spilling.

Allows easy transport of instruments, reagents, and apparatus into the field

Chemical-resistant test kit case is molded from durable, easy-to-clean polypropylene

Technické údaje

Rozmery (V x Š x H):	574 mm x 325 mm x 295 mm
----------------------	--------------------------