



Beaker, griffin, low form, glass, 100 mL

Kat. číslo produktu:

50042H

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Technické údaje

Jednotky:	1 each
Kapacita:	100 mL
Materiál:	Glass, Griffin, Low Form
Materiál skrine:	Glass