



Beaker, griffin, low form, glass, 50 mL, 12/pk

Kat. číslo produktu:	50071H
EUR Cena bez DPH:	Kontaktujte nás
Dátum odoslania neuvedené	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Technické údaje

Jednotky:	12 /pk
Kapacita:	50 mL
Materiál:	Glass, Griffin, Low Form
Materiál skrine:	Glass