



Be Right™



Flask, filtering, 1000 mL

Kat. číslo produktu:

54653

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov

Glass filtering flask with side-arm at 90° angle. Heavy-walled with graduated capacity scale. Designed to ASTM E1406, Type II, Class I requirements. Hose connection accepts 5/16" I.D. flexible tubing.

Technické údaje

Jednotky:	1 each
Kapacita:	1000 mL
Materiál skrine:	Glass