



Recharge kit of 4 pre-filled cartridges with premounted 2956A membranes for A1100 oxygen sensors with 33051-xT cap only

Kat. číslo produktu:

2956A-AT05

Nebezpečné

EUR Cena bez DPH:

Kontaktujte nás

Dodávka do 2 týždňov



Includes o-rings and cleaning tools. Special low concentration electrolyte.

Technické údaje

Parameter:	O2
Podmienky skladovania:	Store in a dry place at room temperature
Prístroj:	Orbisphere A1100-S00T O2 sensor, Orbisphere GA2x00-S00T O2 sensor