



Reagent Kit for EZ1102 Ammonia Analyser

Kat. číslo produktu: APPC1102KTO
EUR Cena bez DPH: Kontaktujte nás
Dátum odoslania neuvedené

Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Technické údaje

Obsah balenia:	APPC1102-01 Ammonia Chlorine Solution (2x)
	APPC1102-02 Ammonia Buffer Solution (2x)
	APPC1102-03 Ammonia Indicator Solution (2x)
	2354153 Ammonia Standard Solution (1x)
Parameter:	Ammonia

Obsah balenia

APPC1102-01 Ammonia Chlorine Solution (2x)APPC1102-02 Ammonia Buffer Solution (2x)APPC1102-03 Ammonia Indicator Solution (2x)2354153 Ammonia Standard Solution (1x)