



## Reagent Kit for EZ1028 Nitrite Analyser (Reagent 1 & 2)

|                      |                 |
|----------------------|-----------------|
| Kat. číslo produktu: | APPC1028KTO     |
| EUR Cena bez DPH:    | Kontaktujte nás |
| Dostupné             |                 |

### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Technické údaje

|                |                                            |
|----------------|--------------------------------------------|
| Obsah balenia: | 52022 SULFANILAMIDE REAGENT 1, 400ML (10x) |
|                | 52023 SULFANILAMIDE REAGENT 2 (10x)        |
| Parameter:     | Nitrite                                    |

---

### Obsah balenia

52022 SULFANILAMIDE REAGENT 1, 400ML (10x) 52023 SULFANILAMIDE REAGENT 2 (10x)