



Test kit, nitrogen, nitrate color disc, model NI-14

Kat. číslo produktu:
EUR Cena bez DPH:
Dátum odoslania neuvedené

1416100-EU
Kontaktujte nás

Nebezpečné



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO₃ and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Model:	NI-14
Model testovacej súpravy:	NI-14**
Najmenší krok:	0,05
Názov metódy:	Color disc/ Cadmium reduction
Parameter:	Nitrogen, Nitrate - As NO ₃ ⁻ -N
Platforma:	Farebný disk
Počet testov:	50
Rozsah merania:	0 - 1 mg/L NO ₃ ⁻ -N
Rozsah merania (2):	0 - 10 mg/L NO ₃ ⁻ -N
Smallest increment:	0,05 0,5
Štandardný rozsah merania:	0 - 1 mg/L NO ₃ ⁻ -N; 0 - 10 mg/L NO ₃ ⁻ -N
Štýl kufrička:	B