



Nitrate test kit, model NI-11

Kat. číslo produktu:

146803-EU

EUR Cena bez DPH:

Kontaktujte nás

Dátum odoslania neuvedené

Nebezpečné



One of Hach's most popular test kits, it uses the simple, low-cost color disc method.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Technické údaje

Model:	NI-11
Model testovacej súpravy:	NI-11**
Názov metódy:	Color disc/ Cadmium reduction
Parameter:	Nitrogen, Nitrate - As NO_3^- -N
Platforma:	Farebný disk
Počet testov:	100
Rozsah merania:	0 - 40 mg/L NO_3^- -N
Smallest increment:	1.0
Štandardný rozsah merania:	0 - 40 mg/L NO_3^- -N
Štýl kufríka:	B